

U.S. Army Research, Development and Engineering Command



Laser-induced plasma chemistry of the explosive RDX with various metals

Jennifer L. Gottfried
U.S. Army Research Laboratory (WMRD)
Aberdeen Proving Ground, MD

TECHNOLOGY DRIVEN. WARFIGHTER FOCUSED.

Presented at: NASLIBS 2011

18 July 2011

Clearwater Beach, FL

UNCLASSIFIED

| maintaining the data needed, and c including suggestions for reducing | lection of information is estimated to ompleting and reviewing the collect this burden, to Washington Headqu uld be aware that notwithstanding an DMB control number. | ion of information. Send comments arters Services, Directorate for Info | s regarding this burden estimate ormation Operations and Reports | or any other aspect of the 1215 Jefferson Davis | nis collection of information, Highway, Suite 1204, Arlington | | |
|---|---|--|---|--|--|--|--|
| 1. REPORT DATE 18 JUL 2011 | | 2. REPORT TYPE | | 3. DATES COVERED 00-00-2011 to 00-00-2011 | | | |
| 4. TITLE AND SUBTITLE | | | 5a. CONTRACT NUMBER | | | | |
| Laser-induced plas metals | 5b. GRANT NUMBER | | | | | | |
| metals | 5c. PROGRAM ELEMENT NUMBER | | | | | | |
| 6. AUTHOR(S) | | 5d. PROJECT NUMBER | | | | | |
| | | 5e. TASK NUMBER | | | | | |
| | | 5f. WORK UNIT NUMBER | | | | | |
| | ZATION NAME(S) AND AE ch Laboratory (WM | 8. PERFORMING ORGANIZATION REPORT NUMBER | | | | | |
| 9. SPONSORING/MONITO | RING AGENCY NAME(S) A | | 10. SPONSOR/MONITOR'S ACRONYM(S) | | | | |
| | | 11. SPONSOR/MONITOR'S REPORT NUMBER(S) | | | | | |
| 12. DISTRIBUTION/AVAIL Approved for publ | LABILITY STATEMENT ic release; distributi | on unlimited | | | | | |
| 13. SUPPLEMENTARY NO | OTES | | | | | | |
| 14. ABSTRACT | | | | | | | |
| 15. SUBJECT TERMS | | | | | | | |
| 16. SECURITY CLASSIFIC | | 17. LIMITATION OF ABSTRACT | 18. NUMBER OF PAGES | 19a. NAME OF RESPONSIBLE PERSON | | | |
| a. REPORT unclassified | b. ABSTRACT unclassified | c. THIS PAGE unclassified | Same as Report (SAR) | 27 | | | |

Report Documentation Page

Form Approved OMB No. 0704-0188



Outline



Motivation

Understand chemistry between metal nanoparticles and molecular explosives

Develop more efficient explosives

LIBS

Nanoparticle production followed by laser-induced plasma chemistry

Time-resolved emission spectra

Laser Parameters

Laser pulse energy dependence

Single vs. double pulse

Substrate characterization

Matrix effects

Effect of impurities on chemistry

Experiments

Double pulse (air/argon)

Aluminum powder additives

RDX discrimination

RDX residue on various metals

Effects of substrate on discrimination

TECHNOLOGY DRIVEN. WARFIGHTER FOCUSED.



Motivation



 IOP PUBLISHING.
 JOURNAL OF PHYSICS D: APPLIED PHYSICS

 J. Phys. D; Appl. Phys. 40 (2007) 3541–3544
 doi:10.1088/0022-3727/40/11/043

Optical investigation of shock-produced chemical products in pseudo-aluminized explosive powders explosion

Ye Song^{1,2}, Wu Jing-he^{1,2}, Wang Yan-ping², Wu Guo-dong² and Yang Xiang-dong^{1,3}

Investigated formation of carbon in aluminized-RDX shock tube initiation

Production of AIO and C decreases rocket propellant performance

XRD confirmed presence of C and Al₂O₃ in blast residue

Observed that ↑ micron-Al results in ↑ C₂ emission

AIO and C₂ emission collected with two 2-nm resolution monochrometers/PMTs

Our idea was to use Al nanoparticles rather than conventional micron-sized Al and look at the emission from additional atomic and molecular species, all on a much smaller scale

¹ Institute of Atomic and Molecular Physics, Sichuan University, Chengdu 610065, People's Republic of China

National Key Laboratory for Shock Wave and Detonation Physics Research, Institute of Fluid Physics, China Academy of Engineering Physics, Mianyang 621900, People's Republic of China




Methodology



- In the past decade, laser ablation has increasingly been used to produce metallic nanoparticles
 - the size and distribution of the nanoparticles can be controlled by:
 - sequential laser pulses, varying repetition rates
 - laser fluence, wavelength and pulse width
 - carrier gas (air, argon, nitrogen, etc.)

J. Laser Micro. Nanoengineer., 3(2), 100-105 (2008).





Methodology



Key advantages of LIBS over the technique used by Song et al.:

Little or no sample preparation is needed

- no need to cast explosive formulations
- any type of material can be ablated with the laser as long as the laser energy > breakdown threshold

The properties of the laser (pulse energy, wavelength, pulse duration) can be tuned to control the size of the particles formed

• laser-generated nanoparticle formation

The intermediate chemical reactions of RDX and Al can be studied on a smaller scale

no shock tube needed

Ability to track relative concentrations of a large number of atomic and molecular species simultaneously


•
$$C + C \rightarrow C_2$$

• C +
$$O_2 \rightarrow CO$$
 + O, CO + O $\rightarrow CO_2$, Al + $CO_2 \rightarrow AlO$ + CO



Laser pulse energy studies








Laser pulse energy studies







Survey of Substrates



Question: do trace impurities in the metal substrate affect the plasma chemistry?

Spectra of 68 metal substrates were acquired

Samples surveyed included:

- high-purity aluminum (99.999%), copper (99.999%), nickel (99.98%), tin (99.998%) and titanium (99.998%)
- numerous metal alloys including brass, lead and steel

Differences in the spectra were observed based on trace element additives and impurities



Aluminum Alloys



Alumina (Al₂O₃)

NIST 1256b (82.99%)


- 9.362% Si (obs.), 3.478% Cu (obs.), 1.011% Zn (obs.)
- <1% Fe, Mn, Ni, Sn, Sr, Ti (obs.)</p>
- <0.1% Cr, Mg, Pb, V (not obs.)

NIST 1259 (89.76%)

- 5.44% Zn (obs.), 2.48% Mg (obs.), 1.60% Cu (obs.), 0.025% Be (obs.)
- <0.2% Cr, Fe, Mn, Ni, Si (not obs.)

NIST 1715 (94.58%)


- 4.474% Mg (obs.), 0.3753% Mn (obs.)
- <0.1% Cr, Cu, Fe, Ni, Pb, Si, Sr, Ti, V, Zn (not obs.)
- Aluminum (99.999%)





Aluminum – C signal








Aluminum – AlO signal







Double pulse experiments



Double-pulse spectra acquired using a Continuum Surelite two-laser system w/echelle spectrograph (EMU-65 with an EMCCD camera)

• 420 mJ per laser, $\Delta t=2\mu s$, $t_{delay}=1.0\mu s$, $t_{gate}=50\mu s$


Spectra of 5 substrates with and without RDX residue were obtained in air and under an argon flow

Al (99.999%), Cu (99.999%), Ni (99.98%), Sn (99.998%), and Ti (99.7%)



RDECOM Double pulse spectra (argon)







Double pulse spectra (argon)




V.I. Babushok, F.C. DeLucia, P.J. Dagdigian, J.L. Gottfried et al., Kinetic modeling study of the laser-induced plasma plume of cyclotrimethylenetrinitramine (RDX), *Spectrochim. Acta, Part B*, **62B**, 1321-1328 (2007).

- the nitrogen must come from the explosive when under argon
- therefore the CN formation is indicative of the chemical reactions the RDX undergoes in the plasma

$$C+N_2 \rightarrow CN+N$$

$$C_2+N_2 \rightarrow 2 CN$$




TECHNOLOGY DRIVEN. WARFIGHTER FOCUSED.



Double pulse results (argon)








RDECOM Double pulse results (argon)









Double pulse results (argon)











<75 µm Al powder mixed with RDX in varying concentrations


- 1:0, 1:1, 1:2, and 1:4 RDX to Al mixtures
- substrates: Al, Cu, Ni, Sn, and Ti
- double pulse laser system under argon

RDX/Al mixtures were crushed onto the substrate surfaces

- each laser shot would blow off a significant amount of material
- 15 spectra of each sample were acquired















- increasing the Al:
 - − ↓ C, H, N, O and CN emission (less RDX sampled?)
 - ↑ C₂ emission (O being scavenged by the Al, less CO and CO₂)







- same optical emission trends observed on all 5 substrates (Al, Cu, Ni, Sn, and Ti)
- ↑ C₂ and AIO confirms observations of Song et al.!
 - decrease in atomic C is new information



RDECOM Classification of RDX on metals



PLS-DA model: 11 classes of pure metals with and without RDX

- RDX+(Al, Cu, Ni, Sn, Ti, Au, Mg, Zn, In, Ag) = 1 class
- Al, Cu, Ni, Sn, Ti, Au, Mg, Zn, In, Ag
- 20 latent variables

Test samples: additional spectra from each sample type

400 spectra in test set

Test samples: spectra acquired on metal alloys not in model

- 160 spectra in aluminum alloy test set
- 280 spectra in other metal alloy test set



RDECOM Classification of RDX on metals




92.5% true positives, 2.5% false positives

| | RDX | Ag | Al | Au | Cu | In | Mg | Ni | Sn | Ti | Zn |
|-------------|-----|----|----|----|----|----|----|----|----|----|----|
| AI+RDX (35) | 34 | 0 | 5 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Cu+RDX (35) | 35 | 0 | 0 | 0 | 1 | 0 | 0 | 0 | 0 | 0 | 0 |
| Ni+RDX (35) | 28 | 0 | 0 | 0 | 0 | 0 | 0 | 33 | 0 | 0 | 0 |
| Sn+RDX (35) | 31 | 0 | 0 | 1 | 0 | 0 | 3 | 0 | 0 | 0 | 0 |
| Ti+RDX (35) | 35 | 0 | 0 | 0 | 0 | 0 | 0 | 4 | 0 | 1 | 0 |
| Au+RDX (5) | 4 | 0 | 0 | 3 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Mg+RDX (5) | 5 | 0 | 0 | 0 | 0 | 0 | 4 | 0 | 0 | 0 | 0 |
| Zn+RDX (5) | 5 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
| Ag+RDX (5) | 5 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| In+RDX (5) | 3 | 0 | 0 | 0 | 0 | 3 | 0 | 0 | 0 | 0 | 0 |
| Ag (5) | 0 | 5 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| AI (35) | 0 | 0 | 35 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Au (5) | 0 | 0 | 0 | 5 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Cu (35) | 0 | 0 | 0 | 0 | 35 | 0 | 0 | 0 | 0 | 0 | 0 |
| In (5) | 2 | 0 | 0 | 0 | 0 | 4 | 0 | 0 | 0 | 0 | 0 |
| Mg (5) | 1 | 0 | 0 | 0 | 0 | 0 | 5 | 0 | 0 | 0 | 0 |
| Ni (35) | 2 | 0 | 0 | 0 | 0 | 0 | 0 | 34 | 0 | 0 | 0 |
| Sn (35) | 0 | 1 | 0 | 0 | 0 | 0 | 0 | 0 | 35 | 0 | 0 |
| Ti (35) | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 35 | 0 |
| Zn (5) | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 5 |



RDECOM Classification of RDX on metals








Classification of RDX on Al alloys (not in model)








Classification of RDX on other alloys (not in model)







Conclusions



Time-resolved, broadband emission of chemical species involved in the reaction of RDX and Al were observed

Confirmed observations of Song et al. using the new experimental methodology

plasma chemistry vs. shock tube detonation

Compared pure metal vs. alloys

 trace metals do affect chemistry, so broadband emission detection is extremely important

Demonstrated that laser pulse energy affects chemistry

related to size of laser ablated particles

Despite differences in the plasma chemistry, RDX residue on different metal substrates can be correctly classified with PLS-DA